

# Bioburden Control, Cleaning and Disinfection



NASA Conference

March 25, 2015

Jim Polarine MA.

Larry Zanko

Science & Solutions for Life

# Agenda

- Bioburden in Cleanrooms
  - Operator Contamination
  - Fungal Spore Contamination
  - Bacterial Spore Contamination
- Cleaning and Disinfection
- *In situ* testing case study




- Tim Sandle
- PDA J Pharm Sci and Tech 2011, 65:392-403
- **A Review of Cleanroom Microflora: Types, Trends, and Patterns**
- Examined isolates from 2000-2009 in U.K.
- Grade A/B and C/D



# Review - Microflora in Cleanrooms (U.K.)

## Grade A and Grade B microflora by group, 2001-2009



# Review - Microflora in Cleanrooms (U.K.)

Genus	A/B (6729)	C/D (2500)
<i>Micrococci</i> (and related)	38%	40%
<i>Staphylococci</i>	21%	11%
<i>Bacillus</i> (and related)	13%	10%
<i>Pseudomonas</i> (and related)	<1%	8%
<i>Corynebacterium</i> (and related)	3%	5%
<i>Rhodococci</i>	<1%	N/A
Fungi	N/A	3%



# Operator contamination

- *Staphylococcus*
- *Propionibacterium acnes*



# Agenda

- ✓ Operator Contamination
- ✓ Fungal Spore Contamination
- Bacterial Spore Contamination



# Fungal Spores

- *Penicillium*
- *Aspergillus*
- *Cladosporium*



*Penicillium*, photos: Ann Larson






# Cleanroom Fungi



Courtesy Dan Klein



# Common sources of Spores

- Items brought into the Cleanroom
  - Bags, Boxes, Intervention Equipment, Pallets, Pallet Jacks, Scrubbers, Cart Wheels, Shoes, Shoe Covers
  - Raw Materials



- ISO-7 Cleanrooms
- Action Levels of 10 and picking up >100
  - Engineering Investigating
  - HVAC
  - Duct Work
  - HEPA Filters
  - Cooling Coils
  - Wall Coverings
  - Airflow Vents




# *Penicillium* Investigation

- Entry and Exit Procedures
- Gowning Procedure
- Cart Wheels
- Construction
  - Further Investigation
  - Use of Sporicides
  - Containers in the Cleanroom
  - Cold room Cleaning Procedures
  - Documentation
  - Assignable Cause



# 70% IPA Efficacy Against Molds


**Fungicidal Activity of 70% Isopropyl Alcohol using Time Kill Method**



# H<sub>2</sub>O<sub>2</sub>/PAA RTU Against Molds - Suspension STERIS

Life Sciences

**Fungicidal Activity of H<sub>2</sub>O<sub>2</sub>/PAA RTU using Time Kill Method**



# Agenda

- **Bioburden in Cleanrooms**
  - Operator Contamination
  - Fungal Spore Contamination
  - Bacterial Spore Contamination
- Bacterial Spore Morphology and Efficacy Testing
- *In situ* testing example



# Bacterial Spores



- *Bacillus cereus* group (7 species\*)
- *Bacillus circulans*
- *Paenibacillus glucanolyticus*

\**B. anthracis*, *B. cereus*, *B. pseudomyoides*, *B. mycoides*, *B. thuringiensis*,  
\**B. weihenstephanensis*, *B. manliponensis*





# Bacterial Endospore




Courtesy Dan Klein



# Exosporium – *B. anthracis*


Cote CK et al. 2011. Microbes and Infection 13(14-15):1146-55.




- ISO-7 and ISO-8 cleanrooms
- Process Vessels
  - Source Locations
    - Cleanroom Shoe Cover
    - Fermentor
    - Process Vessels
  - ✓ The Source was a Raw Material



# Bacillus Testing



# Spore testing results




# Sporicidal Application

- H<sub>2</sub>O<sub>2</sub>/PAA Sporicides
- Cart Wheels
- Items entering the cleanroom



# Gaseous Decontamination



***Complementary approaches which form a highly-effective solution to manage bioburden in critical environments***

# Gaseous Decontamination Methods

DECONTAMINATION METHOD	DELIVERY MEDIUM	PERMISSIBLE EXPOSURE LIMIT*	HUMAN CARCINOGEN	EFFICACY	CYCLE TIME (2500FT <sup>3</sup> )	MATERIAL COMPATIBILITY	REPEATABILITY (VALIDATION)
VHP	Vapor	1.0 ppm	No	Good	< 4 hrs	Good	Good
Hydrogen Peroxide (e.g., fogging, ionization, micro-condensation)	Hybrid	1.0 ppm	No	Good	4-8 hrs	Variable	Moderate
Chlorine Dioxide	Gas	0.1 ppm	No	Good	< 4 hrs	Moderate	Good
Formaldehyde	Gas	0.75 ppm	Yes	Good	> 8 hrs	Good	Good

\* Values represent OSHA permissible exposure limit (PEL) for 8-hour time weighted average (TWA) exposure.

\*\* Visible soils must be properly cleaned before VHP® application.





# Why Use VHP?

- ✓ Efficacy (Broad spectrum sterilant)
- ✓ Consistency & Distribution
  - Reach difficult to access surfaces
  - Passes through HEPA filters
  - Kills airborne and surface microbes
- ✓ Excellent Material Compatibility
  - Electronics
  - Metals and common polymers
- ✓ Speed
  - Minimal labor required
  - Easy to validate (24-hr BI)
- ✓ Green Technology
  - Low toxicity
  - No residues
  - EPA approved



# VHP Process Validation

- Biological Indicators
  - Geobacillus stearothermophilus (Strains 7953 or 12980)
  - 6-log for Biodecontamination / Sterilization
- Environmental Monitoring
- Swabbing



\* *G. stearothermophilus* has been proven to be the most resistant organism to  $VH_2O_2$ .



- Limitations
  - Cannot be used while areas are occupied
  - Higher application cost than liquid chemistries
    - 3<sup>rd</sup> Party Service
    - Capital Equipment and training investment
  - “Lazy” gas needs distribution assistance in large space



# Typical Decon Applications

Anything from a biosafety cabinet to a complete facility...



## Laboratory Equipment

- Microbiological Safety Cabinets
- Transfer Chambers
- Isolators
- Incubators
- Centrifuges

## Complete Facilities

- Clean-rooms
- Corridors
- Laboratories
- Containment facilities
- Changing Areas
- Prep Areas
- Analytical Work areas
- Offices
- Lockers
- Air Showers
- Wash / Toilet areas
- Ductwork
- Filters (HEPA)
- Storage Areas
- Service Areas
- Electrical Cabinets



# When to Decontaminate

- Proactive Basis
  - After shutdown or production change (e.g. Pharma)
  - Before shutdown or equipment service (e.g. BSL-3)
  - Bioburden reduction
- Elimination of Known Contamination (Event Response)
- Commissioning / Decommissioning
  - Facility
  - Equipment



# VHP Process Development for Spaceship Applications

**Authors:** Naresh Rohatgi (NASA JPL) & STERIS Strategic Technology Enterprises

**Publication:** 04ICES-113

**Publication Date:** 2004

- NASA Planetary Protection Office's microbial reduction requirements for all Mars *in situ* life detection missions may require entire spacecraft decontamination.
- Electronics not compatible with approved dry heat methods
- STERIS designed and constructed a high vacuum (~one torr) Biological Indicator Evaluator Resistometer (BIER) vessel to generate hydrogen peroxide lethality data.
- VHP process provided an effective, rapid, safe, and low temperature means for decontaminating spores, mycobacteria, fungi, viruses, and other microorganisms
- VHP process has innocuous residuals as it decomposes to water vapor and oxygen\
- Results under Implementation Plan for Jet Propulsion Laboratory, RG-563852 "Generation of Lethality Data on Vapor Phase Hydrogen Peroxide."



# VHP Certification from NASA Interplanetary Protection

**Authors:** Chen, Fei; Chung, Shirley; Barengoltz, Jack

**Affiliation:** AA(Jet Propulsion Laboratory, California Institute of Technology), AB(JPL), AC(Private Individual)

**Publication:** 38th COSPAR Scientific Assembly, July 2010, in Bremen, Germany, p.4

**Publication Date:** 00/2010

- Flight system must deposit minimal bioload on planets
- Dry heat sterilization (only previously approved method) not suitable for electronics
- Validated using VHP "hardy" strains that were isolated from cleanrooms and environmental populations collected from spacecraft relevant areas.
- Material compatibility discussed



# Case Study: Construction Event at Biotech Site

- ☐ Worst Case Events
- ☐ 9X Clean [1X Sporicide + 2X Phenolic repeated on days 1,2,3]
- ☐ Fogging
- ☐ VHP®
- ☐ Triple Clean
  - ✓ Defined 3X Disinfectants and Sporicide
  - ✓ EM frequency (Static and Dynamic)
  - ✓ Release of the room





# Triple Clean in a Cleanroom

Sample	Action Limit	Pre Triple Clean	Post Triple Clean
RODAC	2 cfu/plate	3 cfu/plate	<1 cfu/plate
RODAC	2 cfu/plate	31 cfu/plate	<1 cfu/plate
RODAC	2 cfu/plate	3 cfu/plate	<1 cfu/plate

Results from an ISO-8 Cleanroom (554ft<sup>2</sup> room)




# In Situ Data-Case Study

Room	Media Type	Action Limits	Pre-Sanitization <sup>a</sup>	Range (#cfu/unit) <sup>b</sup>	Post-Sanitization <sup>a</sup>	Range (#cfu/unit) <sup>b</sup>
#1	Biotest	>2.5 cfu/ft <sup>3</sup>	3 of 4	0.3 <sup>d</sup>	0 of 4	0
	RODAC	>2 cfu/plate	2 of 8	0 to 1	0 of 8	0
	Settling	>2 cfu/plate	0 of 4	0	0 of 4	0
	Swabs	>2 positive	0 of 4	N/A <sup>c</sup>	0 of 4	N/A <sup>c</sup>
#2	Biotest	>2.5 cfu/ft <sup>3</sup>	1 of 4	0.04 <sup>d</sup>	0 of 4	0
	RODAC	>2 cfu/plate	2 of 9	0 to 1	0 of 9	0
	Settling	>2 cfu/plate	0 of 4	0	1 of 4	0 to 1
	Swabs	>2 positive	1 of 7	N/A <sup>c</sup>	0 of 7	N/A <sup>c</sup>




# Cleaning and Disinfection Efficacy




# Time 0

Red = Spore formers


Green = Other




# After 1X Cleaning - NO Sporicide




# After 2X Cleaning – NO Sporicide



# After 3X Cleaning - No Sporicide



# After Sporicide





## Summary

- Bioburden in Cleanrooms
  - Operator Contamination
  - Fungal Spore Contamination
  - Bacterial Spore Contamination
- Bacterial Spore Morphology and Efficacy Testing
- *In situ* testing case study



# Thank You for Your Attendance!



[jim\\_polarine@steris.com](mailto:jim_polarine@steris.com)  
[larry\\_zanko@steris.com](mailto:larry_zanko@steris.com)

